

Ultra**TRAK**TM *ULTIMATE*

REVOLUTIONIZING DIABETIC TESTING

Combines cutting-edge style and performance in an innovative new meter



- Absolutely No Coding
- Strip Ejection Technology
- Results in 5 Seconds and 0.5ul Blood Sample Size
- Brilliant Backlit Display with Glowing Numbers
- Radiant Strip Port Light
- Underfill Detection and 4 Reminder Alarms
- 1000 Test Memory with Date/Time
- Uses AAA Batteries with 5000 test life
- Does *NOT* Use the Dangerous PQQ Enzyme
- 7.14.21.28.60.90 Day Averaging
- Stores Control Solution Testing Results
- Ketone Warning with Pre and Post Meal Record
- Sleek Design, Simple One Button Operation
- Measures Capillary, Venous, and Arterial Blood
- PC Downloadable Results
- Filters Out Most Medications


Product Name	Product Number
Ultra TRAK <i>ULTIMATE</i> Meter	UTUBGM
Ultra TRAK <i>ULTIMATE</i> Test Strips 50/box	UTUBGS
Ultra TRAK <i>ULTIMATE</i> Control Solution (Low and High Levels)	UTUPCS



Manufactured by Vertex Diagnostics. Please visit www.vertexdiagnostics.com for more information.

**Innovating Diabetic Testing So You Can
*Live Your Life***

Live your Life



Ian Brenner, Age 54
Underwater Photographer and Diabetes Patient

**Managing Your Diabetes Just Got Easier,
So You Can Live Your Life to the Fullest**

UltraTRAK™ ULTIMATE

Blood Glucose Monitoring System

Advanced · Sleek · Innovative

No Coding Necessary
For the Simplest Testing & Accuracy

Strip Ejection
For the Safest Handling

Uses AAA Batteries
Easy to Replace

Cutting-Edge Design
The Sleekest Meter on the Market



Rapid Results
In Just 5 Seconds

Tiny 0.5ul Blood Sample
Virtually Painless Experience

Brilliant Backlit Display
Effortless to Read

5000 Test Determinants
Extensive Battery Life

The **UltraTRAK™ ULTIMATE** is the most technologically advanced system on the market. And with competitive pricing, unparalleled customer service, and a five year warranty, we make it easier for you to experience effective diabetes management.




For more information, please visit www.vertexdiagnostics.com or call 888.90.VERTEX

Ultra *TRAK*TM *ULTIMATE*

BLOOD SAMPLE SIZE COMPARISON



RESULTS TIMING COMPARISON















Innovating Diabetic Testing So You Can
Live Your Life

Ultra *TRAK*TM *ULTIMATE*

Delivering the Highest Level of Performance



Competitive Comparisons

	Ultra <i>TRAK</i> <i>ULTIMATE</i> 	Mckesson TrueTrack 	Assure 4 	Assure Pro 	Prodigy Eject 	Medline Optium 
Strip Eject	Yes	No	No	Yes	Yes	No
Autocode	Yes	No	No	No	No	No
Blood Sample Size	0.5ul 	1.0ul 	1.5ul 	1.0ul 	0.6ul 	0.6ul 
Test Speed	5 seconds	10 seconds	10 seconds	10 seconds	6 seconds	5 seconds
Memory	450 tests	365 tests	50 tests	250 tests	450 tests	450 tests
Date/Time	Yes	Yes	Yes	Yes	Yes	Yes
Backlit Display	Yes	No	No	No	No	No
Day Averaging	7.14.30.60.90	14.30	None	7.14.30	7.14.30.60.90	7.14.28
Battery Type	AAA	CR2032	AAA	AAA	CR2032	CR2032
Battery Life	5000 tests	1000 tests	3000 tests	5000 tests	500 tests	1000 tests
Pricing	\$	\$\$\$	\$\$\$	\$\$\$\$	\$	\$\$\$\$

*Date acquired from each manufacturer and up-to-date as of 1/1/11